

Additional Rotors:



Item: Z446-0450
4 x 500ml Angle Rotor (24°)
Max Speed: 8,000 rpm / 10,367 x g
Accel: 59 sec. / Decel: 81 sec.

Inserts/Adapters
Z446-500-A15: 15ml conical
Z446-500-A50R: 50ml round
Z446-500-A50: 50ml conical
Z446-500-A250: 250ml bottle



Item: Z366-06250
6 x 250ml Angle Rotor (28°)
Max Speed: 8,000 rpm / 10,016 x g
Accel: 53 sec. / Decel: 82 sec.

Inserts/Adapters
Z366-06250-A15: 4 x 15ml conical
Z366-06250-A50R: 50ml round
Z366-06250-A50: 50ml conical



Item: Z326-1050
10 x 50ml Angle Rotor (38°)
Max Speed: 10,500 rpm / 16,022 x g
Accel: 56 sec. / Decel: 64 sec.

Inserts/Adapters
Z326-50-A30: 30ml round
Z326-50-A15R: 15ml (conical or round)
Z206-50-A15: 15ml conical
Z206-50-A57: 5-7ml (12-13mm)
Z206-50-A57C: Cushion for 5ml tube



www.hermleUSA.com

PH: 908-769-5555 EM: info@benchmarkscientific.com

Also Available from Benchmark:



Z207 M

High Speed Microcentrifuge

- EZ-Scroll™ touch pad
- High speed, up to 17,317 x g
- Compact, 11in. wide
- Multiple rotor options for 1.5/2.0ml tubes



18 x 1.5/2.0ml
spin column



24 x 1.5/2.0ml
hermetically sealed