

Triopack

new

Acura® manual

- ✓ Three adjustable micropipettes
- ✓ Large choice of volumes
- ✓ Cost-effective prices

Product characteristics

- Exceptional, user-friendly ergonomics
- Ultra-soft, effortless activation
- Permanent volume visibility
- Innovative tip ejection concept
- *Swift-set* user calibration system
- Autoclavable fully assembled at 121°C
- CE certified (IVD 98/79 EEC)

Ordering information

Cat. No.	Contents (Micro / macropipette models)		
825.TRIO.C	2 µL (0.1 - 2)	10 µL (0.5 - 10)	50 µL (5 - 50)
825.TRIO.J	2 µL (0.1 - 2)	20 µL (2 - 20)	200 µL (20 - 200)
825.TRIO.M	10Y µL (1 - 10)	50 µL (5 - 50)	200 µL (20 - 200)
825.TRIO.P	10Y µL (1 - 10)	100 µL (10 - 100)	1000 µL (100 - 1000)
825.TRIO.S	20 µL (2 - 20)	200 µL (20 - 200)	1000 µL (100 - 1000)
825.TRIO.T	50 µL (5 - 50)	200 µL (20 - 200)	1000 µL (100 - 1000)
825.TRIO.V	100 µL (10 - 100)	1000 µL (100 - 1000)	5 mL (0.5 - 5)
835.TRIO.Y	1000 µL (100 - 1000)	5 mL (0.5 - 5)	10 mL (1 - 10)



Each Triopack includes:

- Three pipettes Acura® manual
- Samples of pipette tips
- Operating instruction
- Individual control certificates

Refer to product literature and package insert for more details, and regarding other models.